



# Laser-Induced Damage Threshold (LIDT) Measurement Report

## ISO 21254-2: S-on-1 Test Procedure

Sample: 2-CPW-LO-L/2-1030

**Request from:**

ALTECHNA Co.Ltd.  
Mokslininku st. 6A  
LT-08412 Vilnius

Contact person:

Viktorija Plerpaitė

**Testing institute:**

Lidaris Ltd.  
Saulėtekio al. 10,  
LT-10223, Vilnius, Lithuania, EU

Tester/date:

E. Pupka / 2015-01-14

**Specimen**

Name of sample:

2-CPW-LO-L/2-1030

Type of specimen:

Crystal, AR Coating

Storage, cleaning:

Plastic box, dust blow off by compressed air

**Test specification**

First harmonic of pulsed Nd:YAG InnoLas Laser: SpitLight Hybrid laser ( $\lambda = 1064$  nm, linear polarization, pulse duration 10 ns),  $\lambda/2$  plate combined with additional polarizer attenuator, online scattered light damage detection, offline inspection of damage detection using Nomarski microscopy (100x).

**Laser parameters**

Wavelength:

1064 nm

Angle of incidence:

0 deg.

Polarization state:

linear

Pulse repetition frequency:

100 Hz

Spatial beam profile in target plane:

$TEM_{00}$

Longitudinal beam profile:

Single mode (SLM)

Beam diameter in target plane<sub>(1/e<sup>2</sup>)</sub>:

235.5  $\mu\text{m}$  (average from 64 pulses)

Pulse duration:

10 ns





Fig. 1 Spatial beam profile in target plane (left) and oscilloscope trace (right)

**Test procedure:**

Number of sites per specimen:  
 Arrangement of test sites:  
 Minimum distance between sites:  
 Damage detection:  
 Storage of the specimen:  
 Test environment:  
 Cleaning:  
 Definition of LIDT:

**S-on-1 test**

411  
 Equally spaced  
 744 μm  
 Scattered light diode  
 Plastic box  
 Industrial environment  
 Compressed air  
 Nonlinear fit to 0% of damage probability

Test result:

Table 1 Summarized LIDT's for 2-CPW-LO-L/2-1030

Test mode	Threshold, J/cm <sup>2</sup>
1-on-1	54.50 ≤ 60.59 ≤ 67.01
1000-on-1	21.56 ≤ 24.85 ≤ 27.85





Fig. 2 Characteristic damage curve.

**Typical damage morphology:**



**Fig. 3 Typical in volume damage morphology  
(Energy density 62.39 J/cm<sup>2</sup>, damage after 1 pulse)**



**Fig. 4 Typical front surface damage morphology  
(Energy density 25.98 J/cm<sup>2</sup>, damage after 1000 pulses)**